

	No	Type		No	Type		No	Type
<b>AB Dental</b>			<b>BEGO</b>			<b>Brain Base</b>		
Narrow Platform (NP)	55021	11	Semados Mini	55012	2	Mytis Arrow B-Type 3.7-4.6	55014	4
Internal Hex Platform (SP)	55049	41	Semados S / RI 3.25-5.5	55034	26	Mytis Arrow BS-type 3.7-4.6	55014	4
Conical Platform (CP)	55046	38	Semados RS / RSX 3.0	55082	74	Mytis Arrow C-Type 3.3	55014	4
			Semados RS / RSX 3.75-5.5	55034	26	Mytis Arrow C-Type 3.7-4.6	55020	10
			Semados SC / SCX 3.25-5.5	55034	26	Mytis Arrow E-Type 3.3-4.6	55014	4
						Mytis Arrow Me-Type 4.0-5.0	55020	10
						Mytis Arrow St-Type 4.0-5.0	55020	10
<b>ADIN</b>			<b>Bio3</b>					
Touareg S/OS	55034	26	Advanced and Progressive 3.3-3.8	55025	16			
Touareg Closefit UNP 2.75	55059	51	Advanced and Progressive 4.2-5.0	55067	59			
Touareg Closefit NP 3.0	55059	51				<b>Bredent</b>		
Touareg Closefit RP 3.5	55046	38				SKY	55025	16
Touareg CloseFit WP 4.3, 5.0	55020	10				copaSKY	55027	19
Swell 3.3	55034	26				mini2SKY	55033	25
			<b>BioHorizons</b>			fast & fixed Abutment	55083	75
			Laser-Lok 3.0 Grey	55034	26			
<b>AIRES</b>			Internal Laser-Lok 3.5 Yellow	55049	41			
3.3	55049	41	Internal Laser-Lok 4.5 Green	55051	43			
			Internal Laser-Lok 5.7 Blue	55051	43	<b>BTI</b>		
			External 4.0 Platform	55011	1	Interna Universal	55051	43
						Interna Wide	55051	43
<b>Alfa Gate</b>								
Bioactive, All diam	55034	26	<b>BioMate</b>					
MAX, All diam	55034	26	BioMate 3.3	55074	66	<b>BTK</b>		
M+ NP 3.5	55044	36	BioMate 4.8	55081	73	External Hex EN	55038	30
M+ RP 4.3 - 5.0	55037	29				External Hex ER	55011	1
Porous, All diam	55034	26	<b>Biomet/3i</b>			External Hex EW	55011	1
S-Line 3.0	55059	51	External Hex	55011	1	Internal Hex IR	55049	41
			Certain internal Purple	55036	28	Internal Hex IM	55051	43
<b>Alpha Dent</b>			Certain internal Purple Platform 3.4	55036	28	Internal Hex IW	55051	43
Conical connection Narrow Platform	55025	16	Certain internal Blue, Yellow, Green	55066	58	Hexagonal Morse Taper KR	55081	73
Conical connection Wide Platform	55057	49	T3 Purple	55036	28	Hexagonal Morse Taper KW	55056	48
Hexagon connection Narrow Platform	55039	31	T3 Blue, Yellow, Green	55066	58	Octagonal Morse Taper SR	55014	4
Hexagon connection Regular Platform	55057	49						
			<b>Bioner</b>			<b>C-Tech</b>		
<b>Alpha-Bio</b>			HIKELT 3.8-4.7	55034	26	EL	55067	59
Internal Hex	55034	26	Top DM 3.5-5.0	55046	38			
Conical Hex 3.2	55065	57				<b>Camlog</b>		
Conical Hex 3.75-5.0	55051	43				Screw-Line, Root-Line, Progressive 3.3	55036	28
			<b>Bionika</b>			Screw-Line, Root-Line, Progressive 3.8, 4.3	55031	23
<b>Anker</b>			Cortilog 4.0	55034	26	Screw-Line, Root-Line, Progressive 5.0, 6.0	55032	24
SB III 3.5	55023	14	Cortilog (CCL) 3.3	55065	57	iSy	55065	57
SB III 4.0-5.0	55020	10	Cortilog (CCL) 3.8	55023	14			
			Cortilog (CCL) 5.0	55075	67	<b>Ceralog</b>		
						Hexalobe ceramic	55076	68
<b>Ankylos</b>			<b>Bionnovation</b>					
C/X	55025	16	Biomorse	55023	14	<b>CLC Scientific</b>		
						Conic 3,5-6,0	55081	73
<b>Anthogry</b>			<b>Biotech Dental/Ids</b>					
Axiom Bone Level 3.4 - 5.2	55065	57	Kontakt, Kontakt-S 3.6 - 5.4	55045	37	<b>Conelog</b>		
Axiom Tri-lobe	55065	57	Kontakt-S+ 4.5-5.5	55045	37	Screw-Line, Progressive 3.3	55066	58
Axiom Tissue Level Narrow, Regular	55050	42	Reflect Aspire	55081	73	Screw-Line, Progressive 3.8, 4.3	55066	58
			Reflect Certus	55019	9	Screw-Line, Progressive 5.0	55075	67
<b>Argon</b>			Reflect Tapered	55034	26			
K3pro Rapid 3.0 - 3.5	55026	17	Reflect Recover	55020	10	<b>Cowellmedi</b>		
K3pro Rapid Platf 3mm	55019	9	Reflect Rapid	55020	10	Inno Int Octa	55014	4
K3pro Rapid Platf 3mm Short	55044	36				Inno Sub Hex	55037	29
			<b>Biotem</b>					
<b>Astra Tech</b>			AR 3.7 (Yellow M)	55066	58	<b>Dentack</b>		
Osseospeed EV 3.0	55058	50	AR 4.2-4.6 (Green R)	55037	29	Internal Hex	55056	48
Osseospeed EV 3.6	55046	38	AR 5.1 (Blue W)	55037	29			
Osseospeed EV 4.2	55057	49	BR 3.3 (Yellow M)	55012	2			
Osseospeed EV 4.8	55055	47	BR 3.75-4.0 (Green R)	55014	4	<b>Dental Tech</b>		
Osseospeed EV 5.4	55048	40	BR 5.0-6.0 (Blue W)	55013	3	Implassic FT3 3.25	55039	31
Osseospeed TX 3.0 S	55064	56	IR 3.5 (Yellow M)	55028	20	Implassic FT3 3.75	55049	41
Osseospeed TX 3.5-4.0 S	55046	38	IR 4.1-4.8 (Green R)	55028	20	Implassic FTP	55019	9
Osseospeed TX 4.5-5.0 S	55017	7	IR 4.8-5.5 (Blue W)	55028	20	Implassic HX3	55011	1
						Implassic TR2	55014	4
			<b>BlueSkyBio</b>			Implassic TRW2	55014	4
<b>Avinent</b>			BIO Max NP Platform, All	55046	38			
External, Coral & Ocean 3.5	55038	30	BIO Max RP Platform 4.3-5.0	55037	29	<b>Dentalis</b>		
External, Coral & Ocean 4.1	55011	1	BIO Int Hex 3.5 Platform 3.7-4.1	55049	41	Regular	55051	43
External, Coral 5.1	55013	3	BIO Int Hex 3.5 Platform 4.7	55051	43			
Internal, Coral & Ocean 3.5, 4.1	55051	43	BIO Int Hex 4.5 Platform 4.7-7.0	55051	43	<b>Dentegris</b>		
Conical, Ocean 3.5	55066	58	BIO One-stage Narrow, Regular, Wide	55014	4	3.3	55034	26
Conical, Ocean 4.1	55004	18	BIO One-stage Short Regular & Wide	55014	4	3.75	55057	49
			BIO Quattro Narrow 3.3	55061	53	4.5	55049	41
<b>B&amp;B Implant</b>			BIO Quattro Regular 4.1 - 5.6	55062	54	3.8/4.5	55049	41
Duravit Slim	55021	11	BIO Conus 12, 3.0 Platform 3.0	55064	56	5.5	55020	10
Duravit 3P	55035	27	BIO Conus 12, 3.5/4.0 Platform 3.5-4.0	55046	38			
Duravit Evolution	55035	27	BIO Conus 12, 4.5 x 6mm	55046	38	<b>Dentis</b>		
			BIO Conus 12, 4.5/5.0 Platform 4.5 - 5.0	55017	7	Cleanlant S-Clean 3.7 - 4.8	55020	10
			BIO Conus 12, 5.0 x 6mm	55046	38	S-Clean OneQ Narrow 3.3	55064	56
			BIO Trilobe 3.5	55016	6			
			BIO Trilobe 4.3-6.0	55017	7			

	No	Type		No	Type		No	Type
<b>Dentsply</b>			<b>Hahn</b>			<b>IRES</b>		
PrimeTaper 3.6	55046	38	Tapered 3.0	55058	50	External Hex Narrow 3.3	55012	2
PrimeTaper 4.2	55057	49	Tapered 3.5	55046	38	External Hex Regular 3.7 - 5.2	55011	1
PrimeTaper 4.8	55055	47	Tapered 4.3	55046	38	Internal Hex Narrow 3.3	55087	79
PrimeTaper 5.4	55055	47	Tapered 5.0	55037	29	Internal Hex Regular 3.7 - 5.2	55049	41
			Tapered 7.0	55048	40	Internal Octagon Regular 3.7 - 5.2	55014	4
<b>Dentium</b>			<b>Hi-Tec</b>			<b>J Dental Care</b>		
Superline	55037	29	Tapered Self Thread NTI 3.3, SLT 3,75	55034	26	JDEvolution 3.7-6.0	55067	59
NR Line	55066	58	Tapered Self Thread GTI 4.2, WTI 5.0	55051	43	JDEvolution Plus+ 3.7-6.0	55067	59
Simple Line II GH 1,6	55014	4	Spark TPI	55034	26	JD Evolution S 3.25	55044	36
			Logic Plus	55034	26	JD ICON ultra, S2.75	55064	56
<b>DEVA</b>			Logic LGI 3.0-3.5	55046	38	Conical Abutment	55080	72
Bone Level 3.5	55064	56	Logic LGI 4.3-5.0	55037	29			
Bone Level 4.5	55051	43						
<b>DIO</b>			<b>Hiossen</b>			<b>Keystone</b>		
SM Submerged 3.8	55049	41	HS System	55028	20	<b>I-HEXMR7 7,0-9,0</b>	<b>55051</b>	<b>43</b>
SM Submerged 5.0	55020	10	ET System	55037	29	Prima / Connex	55049	41
UF II Narrow	55064	56	ET 3.2	55047	39	Prima / Connex	55067	59
UF II Regular	55020	10	ET mini	55015	5	Prima / Genesis Narrow	55041	33
			<b>EK 3,5</b>	<b>55046</b>	<b>38</b>	Prima / Genesis Regular	55067	59
			<b>Ek 4,0-5,0</b>	<b>55044</b>	<b>36</b>	Prima / Genesis Wide	55041	33
<b>Ditron Dental</b>			<b>IDI Evolution</b>			Restore Narrow	55079	71
Ultimate 3.75-6.0	55087	79	Standard implant 3.75-4.0	55004	18	Restore Regular	55011	1
MPI 3.3-6.0	55087	79	B1One High	55080	72	Restore Wide	55013	3
			B1One Low	55033	25	TILOBEMAXX 7.0-9.0	55041	33
<b>DSP Biomedical</b>			<b>IDI Implants</b>			<b>Klockner</b>		
External Hexagon (HE) Standard	55014	4	BIO 3.5-5.0	55037	29	VEGA Micro 3.0	55074	66
Hexagonal Internal Standard (HIS)	55034	26	ID Cam ST 4.2	55010	13	VEGA Narrow 3.5	55065	57
Hexagonal Internal Morse (HIM)	55046	38	ID All	55010	13	VEGA Regular 4.0, 4.5	55034	26
Hexagonal Internal Large (HiL)	55037	29	ID3 3.1	55021	11	KL Narrow 3.3	55012	2
Indexed Morse Taper (CMI)	55025	16				KL Regular 3.7, 4.2	55011	1
			<b>Ids /Biotech dental</b>			KL Wide 4.7	55013	3
<b>Easy Implant</b>			Kontakt, Kontakt-S 3.6 - 4.8	55045	37	SK2	55038	30
Ocean	55046	38	Reflect Aspire	55081	73	Essential Cone	55067	59
Lilas	55037	29	Reflect Certus	55019	9	S3M 3.1	55079	71
			Reflect Tapered	55034	26			
<b>Edierre</b>			Reflect Recover	55020	10	<b>Konmet</b>		
3-lobe	55021	11	Reflect Rapid	55020	10	Classic 4.0, 4.8	55014	4
4-lobe	55025	16				Cylindrical 3.5-4.0	55065	57
			<b>Implance</b>			Cylindrical 4.5-5.5	55057	49
<b>Euroteknika</b>			Tissue level	55014	4	Short Cylindrical 5.0-6.0	55078	70
Internal Platform 3.0	55064	56	Bone level Regular	55037	29	Conical 4.4/3.4	55065	57
Internal Platform 3.5	55081	73	Bone level 3.3	55065	57	Conical 4.8/3.8, 5.5/4.5	55057	49
Naturactis/Naturall+ 3.0	55064	56				<b>Kyocera</b>		
Naturactis 3.5-5.0	55081	73	<b>Implant Direct</b>			FINESIA Bone Level 3.2-3.4 (NP)	55058	50
Uneva	55011	1	Legacy 1/2/3/4 3.0/3.2	55012	2	FINESIA Bone Level 3.7-4.2 (RP)	55034	26
			Legacy 1/2/3/4 3.5/3.7 - 5.7/7.0	55051	43	FINESIA Bone Level 4.7-5.2 (WP)	55014	4
<b>Futura implants</b>			RePlant 3.5	55022	12	FINESIA Tissue Level 3.7-4.7 (RP)	55014	4
PRE 4.0	55014	4	RePlant 4.3 - 5.0	55010	13	FINESIA Tissue Level 4.7-5.2 (WP)	55014	4
PRI 4.0	55037	29	RePlant 6.0	55032	24	POI(HAC) EX 3.7 TP/ST	55085	77
			ReActive 3.5 Platform	55022	12	POI(HAC) EX 4.2 TP/ST	55025	16
<b>galimplant</b>			ReActive 4.3 - 5.0 Platform	55010	13	POI(HAC) 2P 3.7/4.2	55004	18
iPx 4012	55046	38	Interactive 3.0 Platform	55046	38			
			Interactive 3.4 Platform	55037	29	<b>Lasak</b>		
<b>GC Implant</b>			Screwindirect 5.0 PF 3.0 - 5.7	55070	62	BioniQ QR	55066	58
Aadva Narrow	55081	73				BioniQ QN	55021	11
Aadva Regular	55066	58	<b>Implant One</b>					
Aadva Wide	55044	36	100 Series Std thread	55029	21	<b>Leader</b>		
Aadva short 4.2	55044	36	200 Series Std thread	55029	21	Internal Hex, 3.3 - 5.5 mm	55034	26
Genesis Plus 3.4-5.0	55085	77	300 Series Std thread (blue)	55029	21	Cono-In	55067	59
SETIO Plus 3.8	55011	1	400 Series Std thread (purple)	55035	27			
			500 series Std thread (green)	55035	27	<b>Maxillent</b>		
<b>Geass</b>			300 Series Wide thread (blue)	55029	21	iRaise	55049	41
Way Extra	55081	73	400 Series Wide thread (purple)	55035	27			
						<b>Medentika</b>		
<b>Global D</b>			<b>Implant Swiss</b>			Microcone NI	55021	11
In-Kone 4.0-5.0	55024	15	Tissue level	55014	4	Microcone RI	55081	73
TwinKon 3.5-4.5	55122	88	Bone level Regular	55037	29	Quattrocone RI	55081	73
			Bone level 3.3	55065	57	Quattrocone 30 AI	55081	73
<b>gmi</b>			<b>Implantvision</b>			Procone 3.3	55119	85
avantgard RP 3.75-4.25	55066	58	Internal Hex 3.7-5.7	55034	26	Procone 3.8	55031	23
avantgard WP 4.75-5.75	55075	67				Procone 4.3	55031	23
frontier RP,WP	55051	43	<b>Innobiosurg</b>			Procone 5.0	55032	24
phoenix NP 3.3	55038	30	MagiCore Regular	55018	8			
phoenix Std 3.3 - Wide 5.0	55011	1				<b>Medentis</b>		
Monolith 3.0	55043	35	<b>Intoss Anchor</b>			ICX Premium/Active Master 3.75 - 4.8 (BL)	55062	54
			Internal Hex 3.7-5.7	55088	80	ICX Premium/Active Master 3.3 (BL)	55058	50
<b>GP Implants</b>						ICX Plus 3.45 Tissue Level	55039	31
SSI Spiral Shape implant	55051	43	<b>Intra-lock</b>					
			Intrahex 3.4 - 4.5	55051	43	<b>Medical Instinct</b>		
<b>gt medical</b>			Gold & Blue 3.4 - 6.0	55051	43	Hex 3.0	55071	63
HXE	55011	1	InDex 3.4 - 3.75	55015	5	Hex 3.4, 4.0, 5.0	55052	44
HXE RP	55012	2	InDex 4.0 - 4.3	55020	10	Cone 3.0, 3.4, 4.0	55071	63
HXE WP	55013	3	FlatTop 3.75	55012	2	Cone 5.0	55072	64
HXI 4.1	55037	29						
002 HXI 4.7	55049	41						
OCT 4.3	55014	4						

	No	Type		No	Type		No	Type
<b>Medigma Biomedical</b>			<b>Nobel Biocare</b>			<b>PLATON JAPAN</b>		
Mars 3,5-6,0	55057	49	Brånemark System Ext hex RP	55011	1	Eight-Lobe Pro NC 3,5-4,1	55046	38
FixTite 3,5-6	55057	49	Brånemark System Ext hex NP	55012	2	Eight-Lobe Pro RC 4,6-5,0	55014	4
Fix-a-Dent 3,75-6	55057	49	Brånemark System Ext hex WP	55013	3	Type IV PRO	55020	10
			Tri-lobe Connection NP	55022	12			
<b>Medimecca</b>			Tri-lobe Connection RP, WP	55010	13	<b>P-I Brånemark</b>		
PT SLA Chaorum	55047	39	Tri-lobe Connection 6.0	55032	24	HEX 3.3	55115	81
			Conical Connection 3.0	55058	50	HEX 3.75 -5.0	55117	83
<b>Megagen</b>			Conical Connection Narrow	55046	38	AMP 3.3	55039	31
Anyridge	55051	43	Conical Connection Regular	55037	29	AMP 3.75 - 4.8	55118	84
Anyone Internal	55037	29	Conical Connection Wide 5.5	55048	40	MT 3.3 - 4.8	55046	38
Mini	55064	56	N1 Narrow	55127	93	Conical Abutment HEX, AMP, MT	55080	72
MUA	55080	72	N1 Regular	55126	92			
Octa	55121	87	MUA NP & RP	55080	72	<b>Prodent Italia</b>		
			On1 NP, RP, WP	55042	34	Twinner & Prime, Conomet, Collar: Orange	55027	19
<b>Meisinger Implants</b>						Twinner & Prime, Conomet, Collar:: Fuchsia	55051	43
OKTAGON Tissue Level RP, WP	55014	4	<b>Noris Medical</b>			Twinner & Prime, Conomet, Collar: Green	55067	59
OKTAGON Bone Level NC	55061	53	Tuff 3.3 - 6.0	55034	26	Twinner & Prime, Conomet, Collar: Yellow	55020	10
OKTAGON Bone Level RC	55062	54				Prime, Conomet, Collar: Blue	55020	10
			<b>Novem</b>					
<b>Microdent</b>			3.4	55039	31	<b>Qualia</b>		
Universal	55011	1	3.8	55016	6	NeoActive	55046	38
Microdent 2.8 - 3.25	55068	60	4.2, 5.0	55057	49			
Microdent 3.3 - 5.0	55011	1				<b>Quest</b>		
Microdent 5.5 - 8.0	55013	3	<b>NucleOss</b>			Hylock NP	55030	22
Genius 3.0	55025	16	Tpure Bone Level	55004	18	Hylock RP	55066	58
Genius 3.5 - 5.0	55051	43	T6 Bonel Level	55074	66			
Ektos 3.2	55023	14				<b>Resista</b>		
Ektos 3.7 - 8.0	55051	43	<b>OCO Biomedical</b>			WIS 3.0-5.0	55049	41
Trylogic 3.5 - 4.2	55051	43	Engage 3.25	55027	19			
Trylogic 5.0	55075	67	Engage 4.0-6.0	55051	43	<b>Rex Implants</b>		
			ERI 3.25	55069	61	External hex TL 1.8	55060	52
<b>MIS</b>			ERI 4.0-5.0	55014	4	Internal hex TL 2.9	55034	26
Conical connection Narrow (C1, V3)	55081	73	TSI 3.25	55068	60			
Conical connection Standard (C1, V3)	55087	79	TSI 4.0-5.0	55011	1	<b>Ritter</b>		
Conical connection Wide (C1)	55041	33				Multi Unit	55033	25
Seven and M4 Narrow	55027	19	<b>OsseoFuse</b>			Spiral SB/LA Narrow PF 2.9	55081	73
Seven and M4 Standard	55049	41	Hexaplus/ HexaPlus-S Sinus 3.75-7.5	55081	73	Spiral SB/LA Standard PF 3.75	55049	41
Seven Wide 5.0	55051	43	ConicalPlus/ ConicalPlus-S Sinus 3.75-5.25	55015	5			
Connect abutment NP, SP, WP	55115	81				<b>RS Implant</b>		
Multi Unit System	55033	25	<b>Osstem</b>			3.5 - 6.0	55020	10
			ETIII 3.5	55015	5			
<b>myplant GmbH</b>			ETIII 4.0-5.0	55037	29	<b>Schütz Dental</b>		
Myplant two 3,5-4,5	55077	69	US Regular	55011	1	Impla cylindrical	55052	44
Myplant Bio 4,0	55077	69	US Mini	55012	2	IMPLA micro Retention Cone 4.2	55052	44
			US Wide	55013	3			
<b>Neo Biotech</b>			US Wide PS	55011	1	<b>SGS Implant</b>		
IS II 4.0	55020	10	SS Regular, Wide	55028	20	P1, P7 Platform	55049	41
IS III 4.0	55020	10	TS Regular, Reg Short, Diam 6 & 7	55037	29	P1D, P5D, P7D	55049	41
Ektos 3.7 - 8.0	55051	43	TS Mini 3.5	55015	5	P7N	55059	51
Trylogic 3.5 - 4.2	55051	43	TS Mini 3.0	55047	39			
Trylogic 5.0	55075	67				<b>SHOFU</b>		
			<b>Osteo Plus</b>			BIOFIX 3.0	55064	56
<b>MIS</b>			Iridium PF 3.8 & 4.2	55079	71	BIOFIX 4.0	55046	38
Conical connection Narrow (C1, V3)	55081	73	She 3.45, 3.75	55011	1	BIOFIX 4.5	55046	38
Conical connection Standard (C1, V3)	55087	79	Shi 3.75, 4.2	55049	41			
Conical connection Wide (C1)	55041	33	Shi 5.0	55051	43	<b>SIC invent</b>		
Seven and M4 Narrow	55027	19				MUA	55120	86
Seven and M4 Standard	55049	41	<b>OT-Medical</b>			SICace 3.4 - 5.0	55052	44
Seven Wide 5.0	55051	43	OT-F2 3.4, 3.8	55053	45	SICmax, SICtapered 3.7 - 5.2	55052	44
Connect abutment NP, SP, WP	55115	81	OT-F2 4.1	55073	65	SICvantage® max, SICvantage tapered 3.0	55026	17
Multi Unit System	55033	25	OT-F2 5.0	55125	91	SICvantage® max, SICvantage tapered 3.7	55116	82
						SICvantage® max, SICvantage tapered 4.2-5.2	55054	46
<b>myplant GmbH</b>			<b>Overmed</b>			<b>Simply Implants</b>		
Myplant two 3,5-4,5	55077	69	Easyline Easydip	55020	10	Internal Hex	55051	43
Myplant Bio 4,0	55077	69	ETLINE	55034	26			
			<b>Oxy Implant</b>			<b>SIN</b>		
<b>Neo Biotech</b>			Internal	55020	10	Strong SW CM	55023	14
IS II 4.0	55020	10	Kone Conical Narrow	55021	11	Strong SW Int Hex	55035	27
IS III 4.0	55020	10	Kone Conical Regular	55051	43	Unitite CM 2.9	55084	76
			PSK Line Regular	55051	43	Unitite CM 3.5	55077	69
<b>Neodent</b>			External Mini	55079	71	Unitite CM Compact 4.0	55085	77
Cone Morse Connection	55025	16	External Regular	55011	1	Tryon 3.25 - 5.0	55011	1
Grand Morse Connection	55066	58	External Maxi	55011	1			
Titamax WS Connection	55056	48						
GM Micro Conical Abutment	55080	72	<b>Paltop</b>					
			Standard 3.75-5.0	55034	26			
<b>Neoss</b>			Narrow	55027	19			
SP (Ø 3.5 implants and wider)	55004	18	PCA Conical	55066	58			
NP (Ø 3.25)	55027	19						
Access Abutment	55060	52	<b>Phibo</b>					
			Aurea EVO NP	55027	19			
			Aurea EVO RP, WP	55051	43			
			TSA S3	55074	66			
			TSA S4, S5	55034	26			
			TSH S3, S4, S5	55011	1			

	No	Type		No	Type		No	Type
<b>Southern Implants</b>			<b>Trinon</b>					
Tri-Nex 4.3-5.0	55010	13	Q2-Implant	55011	1			
Tri-Nex 3.5	55022	12	QK-Implant	55014	4			
External Hex	55011	1	QZA 3.5 Platform	55049	41			
MAX Ext Hex	55011	1	QZA 4.5, 5.7 Platform	55051	43			
Co-Axis Ext Hex 0°, 12°, 24°, 36°	55011	1						
Piccolo Ext Hex	55050	42	<b>URIS</b>					
IT-Connection	55014	4	Omni Narrow 3.0-3.5	55115	81			
Deep Conical 3.0	55064	56	Omni Regular 4.0-6.5	55020	10			
Deep Conical 3.5-4.0	55046	38						
Deep Conical Tapered 5.0	55017	7	<b>Vulkan</b>					
			Internal Hex	55034	26			
			External Hex NP	55012	2			
<b>Spiral Tech</b>			External Hex RP	55011	1			
ESi, Hex connection	55051	43	Conical con MiP	55058	50			
Ultimate, Hex connection	55051	43	Conical con NP	55046	38			
ESi Conical 3.5	55044	36						
ESi Conical 4.3, 5.0	55037	29	<b>Winsix</b>					
			Internal Hex	55020	10			
<b>Straumann</b>								
Tissue level, Regular & Wide Neck	55014	4	<b>Xgate</b>					
Tissue level, Narrow Neck Crossfit (NNC)	55076	68	Conical 3.3 - 5.0	55034	26			
Bone level / Bone Level Tapered Narrow (NC)	55061	53						
Bone level / Bone Level Tapered Regular (RC)	55062	54	<b>Xive</b>					
Bone level / Bone Level Tapered SC 2.9	55064	56	S-Plus 3.0	55030	22			
Bone Level X (BLX) all platform	55027	19	S-Plus 3.4	55061	53			
Screw retained Abutment NC, RC	55080	72	S-Plus 3.8-5.5	55062	54			
PURE Ceramic 4.1	55047	39	TG Plus 3.4-4.5	55018	8			
<b>Surgikor</b>			<b>Zeramex</b>					
Fixation 3.0	55058	50	XT Implant 4.2, 5.5	55086	78			
Fixation 3.5	55066	58						
Fixation 4.3	55037	29	<b>Ziacom</b>					
Solution 3.5, 4.0	55066	58	Conical connection RP	55046	38			
Solution 4.5, 6.0	55037	29	Conical connection WP	55075	67			
Immediate, all	55049	41	External hexagonal connection NP	55079	71			
Versatile, all	55049	41	External hexagonal connection RP	55117	83			
			External hexagonal connection WP	55014	4			
<b>Sweden &amp; Martina</b>			Internal hexagonal connection NP	55027	19			
Syra, Syra SP 3.9-5.0	55011	1	Internal hexagonal connection RP	55057	49			
Outlink <sup>+</sup> 3.75, 4.10, 5.0	55011	1	Internal hexagonal connection WP	55057	49			
Outlink2 3.3	55079	71						
Shelta 3.8-5.0	55025	16	<b>Zimmer</b>					
Premium Straight & SP 3.3-3.8	55025	16	Screw-Vent, Tapered Screw-Vent 3.7-4.1	55049	41			
Kohno Straight & SP 3.3-3.8	55025	16	Screw-Vent, Tapered Screw-Vent 4.7-5.7	55051	43			
Premium Straight & SP 4.25-6.0	55004	18	SwissPlus 3.8	55070	62			
Kohno Straight & SP 4.25-6.0	55004	18	SwissPlus 4.8	55014	4			
Prama 3.8 - 5.0	55079	71	<b>Spline Twist 3.75-5.0</b>	<b>55038</b>	<b>30</b>			
DAT Connection	55035	27	Eztetic NP 3.1	55065	57			
<b>Tag Dental</b>			<b>Zuga</b>					
Axis	55051	43	All	55027	19			
massif	55051	43						
<b>Tatum</b>								
Integrity 3.7 - 5.0	55015	5						
Integrity 6.0 - 8.0	55040	32						
T implant	55063	55						
<b>TBR</b>								
Infinity PF 3.5, PF 4.0	55004	18						
M Implant	55044	36						
<b>Thommen</b>								
SPI 3.5	55021	11						
SPI 4.0	55065	57						
SPI 4.5-6.0	55066	58						
<b>ticare</b>								
Osseus Mini, Standard	55011	1						
Osseus Maxi	55013	3						
Osseus Maxi ps	55011	1						
Inhex Mini	55021	11						
Inhex Standard	55046	38						
Inhex Maxi	55017	7						
Osseus Quattro Standard	55011	1						
Inhex Quattro Standard	55046	38						
<b>TRATE</b>								
ROOTFORM	55065	57						
<b>TRI Dental implants</b>								
Tissue level Octa	55014	4						
Tri-Vent, all	55034	26						
Narrow	55025	16						